

MyCPAP

Presented With Compliments From
www.MyCPAP.com

NEW

Designed to energize.

Zest™

NASAL MASK

Small, quiet, and effortless to fit and use, the Zest™ Nasal Mask offers superb comfort. The innovative design of the Zest™ Nasal Mask combines the Easy-Clip Silicone Seal and Advanced Air Diffuser with proven technology that makes Fisher & Paykel Healthcare masks perform so well.

Simple to both use and seal, the Zest™ Nasal Mask is your perfect solution for an effortless sleep and an invigorated life!


NEW::
Easy-Clip Silicone Seal

One-step removal of Silicone Seal and Foam Cushion for simple, quick and easy use.




Fisher & Paykel
HEALTHCARE

Nasal Delivery – Zest™




The Zest™ Nasal Mask offers an over-the-nose mask solution that fits better, feels lighter and is easier to use. The mask includes more features that make Fisher & Paykel Healthcare masks so popular. The Zest™ Nasal Mask comes as an Out-of-Box Solution fully assembled and ready for use.

UNIQUE MASK FEATURES FOR GREATER FREEDOM, COMFORT AND COMPLIANCE



FlexiFit™ Technology Auto-contouring fit

With Foam Cushion FlexiFit™ Technology automatically adjusts to a wide range of nasal and facial contours with simple one-step ease of fitting and no complicated adjustments.




Glider™ Strap Anti-leak system

The sliding Glider™ Strap provides freedom of movement while maintaining a seal. Mask force remains equal when the head rotates from side to side, minimizing the occurrence of leaks.



NEW:: Easy-Clip Silicone Seal Effortless performance

The Easy-Clip Silicone Seal allows one-step connection or removal for quick and effortless cleaning. The new design allows for perfect seal and frame alignment ensuring optimal mask performance with each use.



NEW:: Advanced Air Diffuser Quiet usage

The maintenance-free Air Diffuser disperses the airflow away from the sleeper's partner, and minimizes exhaust flow noise.



NEW:: Stretchgear™ Headgear Deluxe mask suspension system

A comfortable and breathable elastic headgear, specially designed to maintain a firm seal.



Glider™ Clip One-step removal


Simple one-step removal of headgear. No ongoing headgear adjustments required.



Hidden Strap Release One-step usage

Easy one-step removal of headgear. Enables reassembly without having to make fitting adjustments.

INTERFACE FIT across Normal Population Distribution*



*According to the Fisher & Paykel Healthcare Anthropometric Database.

CONTENTS

- Complete mask with headgear assembled and ready for use
- Single, one-sized Silicone Seal designed to fit most patients
- User-friendly Fitting Guide

REORDER CODE: 400440

FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL FISHER & PAYKEL HEALTHCARE REPRESENTATIVE:

International Fisher & Paykel Healthcare Limited, 15 Maurice Paykel Place, East Tamaki, Auckland 2013, PO Box 14 348, Panmure, Auckland 1741, New Zealand
Tel: +64 9 574 0100 Fax: +64 9 574 0158 Email: info@fphcare.com Web: www.fphcare.com

Australia Tel: +61 3 9879 5022 Fax: +61 3 9879 5232 **Austria** Tel: 0800 29 31 23 Fax: 0800 29 31 22 **Benelux** Tel: +31 40 216 3555 Fax: +31 40 216 3554 **China** Tel: +86 20 3205 3486 Fax: +86 20 3205 2132

France Tel: +33 1 6446 5201 Fax: +33 1 6446 5221 **Germany** Tel: +49 7182 93777 0 Fax: +49 7182 93777 99 **India** Tel: +91 80 4123 6041 Fax: +91 80 4123 6044 **Irish Republic** Tel: 1800 409 011

Italy Tel: +39 06 7839 2939 Fax: +39 06 7814 7709 **Japan** Tel: +81 3 3661 7205 Fax: +81 3 3661 7206 **Middle East** Tel: +9 71 4 360 0849 Fax: +9 71 4 368 8110 **Northern Ireland** Tel: 0800 132 189

Spain Tel: +34 902 013 346 Fax: +34 902 013 379 **Sweden** Tel: +46 8 564 76 680 Fax: +46 8 36 63 10 **Switzerland** Tel: 0800 83 47 63 Fax: 0800 83 47 54 **UK (EU Authorised Representative)**

Fisher & Paykel Healthcare Limited, Unit 16, Cordwallis Park, Clivemont Road, Maidenhead, Berkshire SL6 7BU, **United Kingdom** Tel: +44 1628 626 136 Fax: +44 1628 626 146

USA Tel: 1800 446 3908 +1 949 453 4000 Fax: +1 949 453 4001